

Introducing



Combination with Style— The Sonicator® Plus 940

Introducing the new Sonicator® Plus 940 four channel combination unit. It combines four channels of electrical stimulation with 1 and 3 MHz therapeutic ultrasound.

There are eight treatment waveforms that include: Interferential, Premodulated, Medium Frequency (Russian), EMS, High Volt, TENS, Microcurrent and Direct Current (DC). All waveforms are available on all four channels.

To get you started quickly, the Sonicator Plus 940 has 61 preset protocols. In addition, there is space to store up to 80 user-defined protocols so that you can customize treatment programs for your practice.

The membrane panel and back-lit LCD touch screen allow you to quickly setup

a treatment. A large variety of options for both ultrasound and electrical stimulation can be chosen to customize the Sonicator Plus 940 for your patients' conditions.

The Sonicator Plus 940 comes with a dual-frequency 5.5 cm² applicator. An optional 0.9 cm² applicator is also available. Both applicators feature a contact indicator that alerts you when coupling is inadequate.

Mettler Electronics celebrates its Fiftieth Anniversary this year, as a family owned and operated company dedicated to provide quality products to you, its customer. To receive a free color brochure describing the Sonicator Plus 940, contact us at 1.800.854.9305

