



FOUR CHANNELS AND ALL THAT JAZZ.

The New Mettler Trio.™

Introducing three instruments that harmonize perfectly with your therapy needs. These flexible, full-featured Mettler systems set a new beat for value—they're so complete, there's really no other 'jazz' to buy. Consider their cool capabilities:

- The four-channel Sonicator® Plus 994 gives you four channels of muscle stimulation plus ultrasound.
- The Sonicator® Plus 992 has two channels of muscle stimulation plus ultrasound.

- The Sys*Stim® 294 provides four channels of muscle stimulation.


Each of Mettler's versatile performers can treat a single or multiple patients with amazing ease. Membrane control panels, soft-touch knobs and LEDs take you smoothly through each therapy session. It's quality you can feel—backed by ISO 9001 Certification and a two-year limited* warranty to protect your investment.

So call your Mettler Electronics representative or (800) 854-9305 today.



METTTLER **ELECTRONICS**®

THE METTLER TRIO™



The Mettler Trio includes the Sonicator Plus 994, the Sonicator Plus 992 and the Sys*Stim 294.



The Sonicator Plus units include a dual frequency, 5cm² applicator. Optional 1cm²/3.3 MHz and 10cm²/1 MHz are also available.

	<i>Sonicator Plus 994</i>	<i>Sonicator Plus 992</i>	<i>Sys*Stim 294</i>
90-240 V, 50/60 Hz input; can be used just about anywhere	✓	✓	✓
Independent timers for different treatment setups	5	3	4
Stimulation channels	4	2	4
Six treatment waveforms; Interferential, Premodulated Medium frequency, Biphasic (TNS), High volt & Microcurrent	✓	✓	✓
Vector rotation, surge and reciprocation	✓	✓	✓
1 & 3 MHz ultrasound	✓	✓	
Interchangeable applicators with universal applicator cable	✓	✓	
100%, 50% and 20% Duty cycles	✓	✓	
Combination therapy	✓	✓	
ETL and C-ETL Listed (9801427)	✓	✓	✓
Year 2000 Compliant	✓	✓	✓



Independent timers allow simultaneous treatment setups for multiple patients or treatment areas.



Add an optional treatment cart to create a mobile treatment center in your office.