



Technical Specifications

The G5 GBM hand-held massager provides 15 to 55 CPS. It is ergonomically designed to reduce hand vibration. Use the G5 GBM massager for home use or light-duty clinic only.

G5 massage machines have been sold and used successfully for over 50 years. G5 massage machines are the only massagers with Directional Stroking[®], not vibration.

The G5 GBM provides very deep penetrating massage. It is reliable, durable and quiet. Using the G5 GBM massage machine over a hot pack or cold pack seems to drive in the heat or cold. It can be equally-effective when used over towelings or clothing, and won't tangle in hair.

Accessories include: 1- four-ball applicator (6009); 1-Two-Ball Applicator (6012); 1- pointed-tip applicator (6014); 1 large 5" diameter firm rubber applicator (6015); 1- curved sponge applicator (6016); 25-medium disposable applicator covers (6038) and; 25-disposable applicator covers (6039).



Motor:	1/40 horsepower
Output Speed:	15 to 55 cycles per second
Power Source:	90 - 240 VAC 50/60 H
Speed Control:	Solid State Electronic
Leakage Current:	under 75 μ A
Unit Dimensions:	10" x 4" x 4"
Unit Weight:	3 Lbs. 10 oz.
Applicator Packages:	PRO-PORT-PKG (60002)
Shipping Weight:	10 lbs.



Mettler Electronics Corp.