

# G5® Applicators & Applicator Covers

Part No.	Description	Where Used
AP209	Right-Angle Directional-Stroking® Adaptor, USA Brand	Deep directional stroking on large muscles
AP209F	French Right-Angle Directional-Stroking® Adaptor (Adjustable)	Variable angle directional stroking action
AP210	Soft Sponge Rubber Applicator, 2 5/8" Diameter	Gentle massage for sensitive body areas
AP212	Soft Sponge Rubber 3 1/2" Diameter	Gentle massage for sensitive body areas
AP214	Applicator with one fine sponge	Gentle surface massage
AP214FC	Coarse Foam Replacement Pads for AP214, Package of 10 pads	
AP21FF	Fine Foam Replacement Pads for AP214, Package of 10 pads	
AP215	Multiple Prong Applicator	Vigorous massage over large body areas
AP216	Four-Ball Firm Rubber Applicator	Deep massage of major muscle groups
AP217	Half-Ball Firm Rubber Applicator	Soles of the feet or hands to relax muscles
AP222	Blunt-Tipped Firm Rubber Applicator	Deep-working of soft tissue
AP223	Two-Ball Firm Rubber Applicator	Used to straddle the spinal column
AP225	Scalp and Skin Surface Applicator	Used when hyperthermia and warmth is desired
AP227	Pointed-Tip Firm Rubber Applicator	Used to reduce trigger points
AP229	Large Firm Rubber Applicator	Deep massage on the back, buttocks and thighs
AP230	Curved Sponge Applicator	Massage of body limbs and the lower back
AP231	Flex Sponge Rubber Applicator Insert for AP230	
AP232	Flexible Applicator Frame for AP230	
AP331	Adapter Ring, To Protect Applicator Threads & Allow Use of all Applicators	
AP400	Head Holding Handles, Control Guidance System for the Macro Head	
AP509	Micro Right Angle Adaptor, USA Brand	Used for small, delicate body areas
AP509F	Micro Right Angle Adaptor French Brand (Adjustable)	Multiple angles for small, delicate body areas
AP510	Micro Round Soft Sponge Rubber	Used on face and other delicate body areas
AP514	Micro Round Fine Sponge with Back (1)	Delicate body areas
AP514FS	Package of 10 Micro Fine Sponge replacements	
AP514CS	Package of 10 Micro Coarse Sponge replacements	
AP522	Micro Blunt-Tipped Firm Rubber	Delicate body areas
AP525	Micro Skin & Scalp Applicator	Scalp and skin stimulation
AP236	6 Plastic Reusable Covers for AP230	
AP234D	50 Small Disposable Covers for AP 210, 217, 250, 251, 252 & 253	
AP235D	50 Medium Disposable Covers for AP 212, 216 & 229	
AP236D	50 Large Disposable Covers for AP230	



1333 South Claudina Street \* Anaheim, CA 92805 \* U.S.A.  
 Toll Free: (800) 854-9305 \* Tel: 1(714) 533-2221 \* Fax: 1(714) 635-7539  
 Web site: [www.mettlerelectronics.com](http://www.mettlerelectronics.com) \* E-mail: [mail@mettlerelectronics.com](mailto:mail@mettlerelectronics.com)