



# Pro Power<sup>®</sup>

## Technical Specifications



Motor -	1/40 horsepower
Output Speed -	20-50 cycles per second
Power Source -	90 - 240 VAC 50/60 Hz
Speed Control -	Solid State Electronic
Current Leakage -	under 75 MicroAmps
Power Cord -	10' Standard 3-Prong Plug
Housing Material -	ABS Cyclolac
Configurations -	Hand-held
Unit Dimensions -	3" H x 12" W x 3" D
Unit Weight -	3 Lbs.
Shipping Dimensions -	17" x 11" x 9"
Shipping Weight -	7 Lbs.
Applicator Packages -	PRO-PORT-PKG



**Mettler Electronics Corp.**