

# 6D ACTION<sup>™</sup>

Six-Dimensional Method

## 6D Action Device

**6D Action Device is based on negative pressure or vacuum technology and can be used as an aid for physiotherapy and lymphatic drainage therapy.**

6D Action vacuum based device allows continuous vacuum suction and also pulsating vacuum suction. The device package includes the main unit, 6 plastic treatment cups in different sizes, single and dual vacuum hoses with a filter for plastic treatment cups, a roller treatment head with a filter, a power cable and a carrying bag.



**6D Action device package**

1. Main unit
2. Single hose for plastic treatment cups + filter
3. Dual hose for plastic treatment cups + filter
4. Plastic treatment cups x 6, diameters 0.98 in, 1.18 in, 1.38 in, 1.77 in, 2.56 in, 3.35 in
5. Roller treatment head + filter, diameter 2.95 in
6. Power cable
7. Carrying bag

### Accessories

8. Glass treatment cups x 5, diameters 2.52 in, 1.85 in, 1.42 in, 0.47 in, 0.43 x 0.75 in (oval)
9. Single hose for plastic treatment cups + filter
10. Dual hose for plastic treatment cups + filter
11. Single hose for glass treatment cups + filter
12. Dual hose for glass treatment cups + filter
13. Filter for single and dual hose

14. Roller suction cup set of 5 filters

15. Roller treatment heads + filter in two sizes, diameter 2.95 in and 1.97 in

### Technical data

Vacuum: 0-300 mm Hg(± 10-15%)  
 Pulsation: 0.1-3.0 s (± 10-15%)  
 Dimensions: 7.1 x 10.2 x 7.5 in  
 Weight: 10.14 lbs

## New!

Autoclavable glass treatment cups and a dual hose for using two treatment cups simultaneously.

Glass treatment cups, dual and single hoses for both plastic and glass treatment cups are now available as individual accessories.



## Six-Dimensional Treatment with 6D Action

Power requirements: 110 VAC, 60 Hz, 90 VA  
 CE, EMC: 2014/30/EU and LVD: 2014/35/EU

[www.6dtape.com](http://www.6dtape.com)



1333 S. Claudina St.  
 Anaheim, CA 92805  
 Toll Free: 800.854.9305  
 Tel: 1.714.533.2221  
[www.mettlerelectronics.com](http://www.mettlerelectronics.com)  
[mail@mettlerelectronics.com](mailto:mail@mettlerelectronics.com)