

APT[®] Specifications

APT Unit

Model 8300

Weight	22 lbs.
Length	28 in. Overall
Width	18 in.
Height	6.5 in. Folded
Working voltage	24 VDC
Revolutions per minute	20 - 60 RPM
Rotation Radius	1.5, 3.0, 4.5, 6.0 in.
Exercising Positions	15°, 30°, 45°, 60°
Current Rating	4.16 A Max.
Accuracy of measurement display	±10%



Type B equipment

Equipment not suitable for use in the presence of flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

Continuous operation

External Power Supply

Weight	1.8 lbs.
Length	7.5 in.
Width	3.5 in.
Height	1.8 in.
Input:	100~240 V 47-63 Hz 1.25 A
Output:	24 VDC 100 W Max. 4.16 A



Class I equipment

APT Hi Lo

Model 8500

Weight	86 lbs.
Base Dimensions	40 x 30 in.
Height	40 in.
Working Height from Ground Up	50 in.



Mettler Electronics Corp.