

CPM 8000

Mettler Continuous Passive Motion

Introducing the Mettler CPM 8000. It provides complete Automatic motion of the knee. It is ideal for all kinds of knee CPM patients including orthopedic, pediatric and bariatric.



SPECIFICATIONS

Input:	100 to 230 VAC, 50/60 Hz
Weight:	33 pounds, (15 kg)
Dimensions:	47" x 12" x 16.5", (120 x 30 x 42 (cm))
Treatment Time:	Continuous
Speeds	
Min. (100%)	70 s \pm 10% (- 10° to 120°) 75 s \pm 10% (120° to - 10°)
Mean (60%)	120 s \pm 10% (- 10° to 120°) 120 s \pm 10% (120° to - 10°)
Max. (30%)	240 s \pm 10% (- 10° to 120°) 250 s \pm 10% (120° to - 10°)
Working Angles:	- 10° (\pm 2°) to 120° (\pm 2°)
Range of Motion: (ROM)	For H = 8 from -10° to 120° For H = 5 from -10° to 115° For H = 0 from -10° to 110°

Features

- Easy to use
- Adjustable
- Patient lock-out
- Force reversal feature
- "Whisper Quiet" screw drive
- Wide carriage for all sizes of patients
- 2 Year Warranty
- Programmable
- Reversible footplate
- Easy to use LCD control unit



Mettler Electronics Corp.

For more information go to
www.MettlerElectronics.com
or call 1 (800) 854-9305

Copyright © 2012 by Mettler Electronics Corp.—Anaheim, CA