



G5 Plus™

Technical Specifications

The G5 Plus is the newest G5 massager in a long line of clinical models. Its operating speed is adjustable from 10 to 60 CPS (cycles per second). The G5 Plus features a speed control dial. The back-lit display shows massaging speed as well as treatment time remaining. It maintains its speed (CPS), regardless of load. Perfect for busy clinics where therapy is often done from both sitting and standing positions. It is especially useful for over-the-shoulder applications.

G5 massage machines have been sold and used successfully for over 50 years. G5 massage machines are the only massagers with Directional Stroking®, rather than “thumping”. The G5 Plus provides very deep, penetrating massage. It is reliable, durable and quiet. Using the G5 Plus massage machine over a hot pack or cold pack seems to drive in the heat or cold. It can be equally-effective when used over towelings or clothing, and won't tangle in hair.

Accessories include: 1-large sponge 3½" diameter applicator (6004); 1- multiple-prong applicator (6008); 1- four-ball applicator (6009); 1- pointed-tip applicator (6014); 1 large 5" diameter firm rubber applicator (6015); 1- curved sponge applicator (6016); 50-medium disposable applicator covers (6038) and; 50-large disposable applicator covers (6039).



Motor:	1/17 horsepower
Output Speed:	10 to 60 cycles per second
Power Source:	90 - 240 VAC, 50/60Hz
Product Approval:	CE
Leakage Current:	under 75 µA
Size:	44" x 21" x 21"
Weight:	25 pounds
Applicator Packages:	PRO -PKG (60001)
Shipping Weight:	35 lbs.
Shipping Dimensions:	22" x 21" x 17.4"



Mettler Electronics Corp.