



Mettler Electronics Corp.



**The Ultimate in
Professional
Massage**



G5® Plus®

The G5 Plus is the newest massager in a long line of clinical models. It features a digital display that shows speed and time. Speed is adjustable from 10 to 60 cycles per second (CPS). The G5 Plus maintains its velocity, regardless of load. "Directional Stroking™" provides very deep penetration that can more effectively deliver hot or cold pack treatments. The G5Plus is ideal for physical therapy clinics and sports medicine. It helps in trigger point reduction, treatment of muscle spasm and pain relief. It's a powerful and quiet machine, perfect for busy clinics.



Item#:	ME610
Motor:	1/17 horsepower
Input:	90 - 240 VAC, 50/60Hz
Assembled Dimensions:	44" x 21" x 21"
Assembled Unit Weight:	25 Lbs.
Shipping Dimensions:	22" x 21" x 17.4"
Shipping Weight:	35 Lbs.
Applicator Package:	PRO-PKG
Approvals:	CE
Warranty:	3 years



G5® GBM®

This hand-held massager provides 15 to 55 CPS. It is ergonomically designed to reduce hand vibration, which is an important feature that allows multiple treatments in the same day without fatigue. The G5 GBM incorporates power and mobility in a small unit that fits in a travel bag. Simple to operate—just turn the knob and gradually go from lower to higher speeds. It is perfect for light-duty clinical treatments or home use.



Item #:	ME600
Motor:	1/40 horsepower, 24 volt DC
Input:	90 to 240 VAC, 50/60 Hz
Assembled Dimensions:	10" x 4" x 4"
Assembled Unit Weight :	3 Lbs.
Shipping Dimensions:	18" x 12" x 9"
Shipping Weight:	10 Lbs.
Applicator Package:	PRO-PORT PKG
Approvals:	UL, TUV, CE, RoH5
Warranty:	1 year



Optional G5 Travel Bag



G5® Therassist®

The G5® Therassist is our latest advancement in massage/percussion technology. This unit features digital solid-state electronic controls that shows the timer and speed control with ranges from 15 to 60 CPS. It more closely replicates manual procedures. Ideal for stress reduction, myofascial release, sports rehabilitation, postural drainage. It comes complete with caster-stand and tray. It can be used as table-top or wall mounted.

Item #:	TRS 24 W
Motor:	1/17 horsepower, 24 volt DC
Input:	90 to 240 VAC, 50/60 Hz
Assembled Dimensions:	44" x 21" x 21"
Assembled Unit Weight:	25 Lbs.
Shipping Dimensions:	22" x 21" x 17 ³ / ₈ "
Shipping Weight:	35 Lbs.
Applicator Package:	PRO-PKG
Approvals:	UL, TUV, CE, RoH5
Warranty:	1 year

Wall bracket and mounting hardware are sold separately



G5® Gemini®

The most versatile model providing physical therapy, respiratory therapy, and massage therapy applications. The speed ranges from 15 to 60 CPS. This model offers a macro head and a micro head for specialized therapy. The micro head is used for pediatrics, geriatrics, and facial areas. The micro head provides concentrated therapy while the macro head provides powerful therapy to larger muscle groups. It comes complete with custom cart with built-in bins for storage and twelve applicators.

Item #:	GEM 24 W
Motor #1:	1/17 horsepower, 24 volt DC
Motor #2:	1/40 horsepower, 24 volt DC
Input:	90 to 240 VAC, 50/60 Hz
Assembled Dimensions:	50 ¹ / ₄ " x 18" x 22 ¹ / ₂ "
Assembled Unit Weight:	48 Lbs.
Shipping Dimensions:	2 boxes, 30" x 23" x 19" and 43" x 19" x 18"
Shipping Weight, 2 boxes:	64 Lbs.
Applicator Packages:	PRO-PKG and MCT-PKG
Approvals:	UL, TUV, CE, RoH5
Warranty:	1 year





G5® GK3®

High quality design for heavy and frequent use. Versatile and easy to start a treatment. Quiet, powerful Directional Stroking® percussor apparatus designed specifically as an adjunct for massage. Basic control panel with ranges from 15 to 60 CPS. It comes with caster-stand and tray, and also can be used as table top.

Item #:	GK 24 B
Motor:	1/17 horsepower, 24 volt DC
Input:	90 to 240 VAC, 50/60 Hz
Assembled Dimensions:	44" x 21" x 21"
Assembled Unit Weight:	25 Lbs.
Shipping Dimensions:	22" x 21" x 18"
Shipping Weight:	36 Lbs.
Applicator Package:	PRO-PKG
Approvals:	UL, TUV, CE, RoH5
Warranty:	1 year



G5® Workout Masseur®

This unit features a continuously-variable speed range of 15 to 60 CPS and control dial with CPS indicator. Ideal for high volume physical therapy settings. It comes complete with mobile caster stand and storage tray. It can be used as pole, table top or wall mounted.

Item #:	WM24W
Motor:	1/17 horsepower, 24 volt DC
Input:	90 to 240 VAC, 50/60 Hz
Assembled Dimensions:	44" x 21" x 21"
Assembled Unit Weight:	25 Lbs.
Shipping Dimensions:	22" x 21" x 18"
Shipping Weight:	36 Lbs.
Applicator Package:	PRO-PKG
Approvals:	UL, TUV, CE, RoH5
Warranty:	1 year

Wall bracket and mounting hardware are sold separately.



G5[®] Pro Power[®]

This massager is beneficial to anyone who works out, does strenuous work or competes. Before an activity, use it to enhance flexibility, increase range-of-motion, and reduce the risk of muscle spasms, strains or contusions. After exertion, rely on the G5 Pro Power massager to accelerate the healing process and to reduce muscle pain caused by work and stress. Variable speed from 15 to 55 CPS, an easy on/off control and a directional applicator, make this unit ideal for self-massage at home.

Item #:	PP24B
Motor:	1/40 horsepower, 24 volt DC
Input:	90 to 240 VAC, 50/60 Hz
Assembled Dimensions:	3" x 12" x 3"
Assembled Unit Weight:	3 Lbs.
Shipping Dimensions:	17" x 11" x 9"
Shipping Weight:	7 Lbs.
Applicator Package:	PRO-PORT-PKG
Approvals:	UL, TUV, CE, RoHS
Warranty:	1 year



G5[®] Applicator Packages

Each G5 massage machine comes with an applicator package. Below are listed the packages and their contents.



PRO-PKG

- Large Sponge Applicator, 3½" Diameter
- Multiple-Prong Applicator
- Four-Ball Applicator
- Pointed-Tip Applicator
- Large Round Firm Rubber, 5" diameter
- Curved Sponge Applicator
- (50) Medium Disposable Applicator Covers
- (50) Large Disposable Applicator Covers



PRO-PORT-PKG

- Four-Ball Applicator
- Two-Ball Applicator
- Pointed-Tip Applicator
- Large Round Firm Rubber 5" diameter
- Curved Sponge Applicator
- (25) Medium Disposable Applicator Covers
- (25) Large Size Disposable Applicator Covers



MCT-PKG

- Micro Right Angle Adaptor, USA Brand
- Micro Soft Sponge Rubber
- Micro Fine Sponge with Back
- (10) Micro Fine Sponge Replacements
- (10) Micro Coarse Sponge Replace.
- Micro Blunt-Tipped Applicator
- Micro Skin & Scalp Applicator
- Micro Round Firm Rubber Applicator
- (50) Extra Small Disposable Applicator Covers

G5® Applicators & Applicator Covers

Part No.	Description	Where Used
AP209	Right-Angle Directional-Stroking® Adaptor, USA Brand	Deep directional stroking on large muscles
AP209F	French Right-Angle Directional-Stroking® Adaptor (Adjustable)	Variable angle directional stroking action
AP210	Soft Sponge Rubber Applicator, 2½" Diameter	Gentle massage for sensitive body areas
AP212	Soft Sponge Rubber 3½" Diameter	Gentle massage for sensitive body areas
AP214	Applicator with one fine sponge	Gentle surface massage
AP214FC	Coarse Foam Replacement Pads for AP214, Package of 10 pads	
AP21FF	Fine Foam Replacement Pads for AP214, Package of 10 pads	
AP215	Multiple Prong Applicator	Vigorous massage over large body areas
AP216	Four-Ball Firm Rubber Applicator	Deep massage of major muscle groups
AP217	Half-Ball Firm Rubber Applicator	Soles of the feet or hands to relax muscles
AP222	Blunt-Tipped Firm Rubber Applicator	Deep-working of soft tissue
AP223	Two-Ball Firm Rubber Applicator	Used to straddle the spinal column
AP225	Scalp and Skin Surface Applicator	Used when hyperthermia and warmth is desired
AP227	Pointed-Tip Firm Rubber Applicator	Used to reduce trigger points
AP229	Large Firm Rubber Applicator	Deep massage on the back, buttocks and thighs
AP230	Curved Sponge Applicator	Massage of body limbs and the lower back
AP231	Flex Sponge Rubber Applicator Insert for AP230	
AP232	Flexible Applicator Frame for AP230	
AP331	Adapter Ring, To Protect Applicator Threads & Allow Use of all Applicators	
AP400	Head Holding Handles, Control Guidance System for the Macro Head	
AP509	Micro Right Angle Adaptor, USA Brand	Used for small, delicate body areas
AP509F	Micro Right Angle Adaptor French Brand (Adjustable)	Multiple angles for small, delicate body areas
AP510	Micro Round Soft Sponge Rubber	Used on face and other delicate body areas
AP514	Micro Round Fine Sponge with Back (1)	Delicate body areas
AP514FS	Package of 10 Micro Fine Sponge replacements	
AP514CS	Package of 10 Micro Coarse Sponge replacements	
AP522	Micro Blunt-Tipped Firm Rubber	Delicate body areas
AP525	Micro Skin & Scalp Applicator	Scalp and skin stimulation
AP236	6 Plastic Reusable Covers for AP230	
AP234D	50 Small Disposable Covers for AP 210, 217, 250, 251, 252 & 253	
AP235D	50 Medium Disposable Covers for AP 212, 216 & 229	
AP236D	50 Large Disposable Covers for AP230	



1333 South Claudina Street * Anaheim, CA 92805 * U.S.A.
 Toll Free: (800) 854-9305 * Tel: 1(714) 533-2221 * Fax: 1(714) 635-7539
 Web site: www.mettlerelectronics.com * E-mail: mail@mettlerelectronics.com